For Research Use Only. Not for use in diagnostic procedures.

HA-tagged Protein Purification Gel

(MoAb. clone 5D8)

Code No. 3321

*Purification to Maintain Protein Activity
from eukaryote cell lysate and culture supernatant*

Product Description

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more. Thus, choosing a method of purification is an important aspect in maintaining protein structure and function.

Recombinant tagged protein purification methods and kits are now widely recognized. The HA epitope tag, YPYDVPDYA, is derived from Human influenza hemagglutinin (HA) 98-106 aa. This HA-tag is also commonly used in such as mammalian cell expression vectors.

MBL's HA-tagged Protein Purification Gel is designed for the isolation of HA-tagged protein from cell culture supernatants and cell lysate under neutral pH condition. Severe conditions such as acidic or alkaline elution denature protein structure. However, a neutral pH elution can preserve protein activity and native conformation. Therefore, MBL has developed the Anti-HA-tag Gel to purify HA-tagged proteins quickly and efficiently at neutral pH to maintain the activity and conformation of purified proteins. The elution of HA-tagged proteins from the Gel is achieved by the addition of the HA-tag peptide (YPYDVPDYA). As the HA-tag peptide competes with HA-tagged proteins on the Gel, the purified proteins do not lose the protein activity.

Components

	Quantity
	Code No. 3321
<u>Anti-HA-tag Gel</u>	1 mL × 1 vials
50% slurry: 1 mL Gel in 2 mL total volume in PBS with 0.09 % sodium azide as preservative.	

Elution Peptide (YPYDVPDYA) is NOT included.

It is available as Code No. 3320-205 (2 mg x 5 vials).

Product Capacity

The purification capacity of the Anti-HA-tag Gel varies depending upon the HA-tagged protein. For example, 1 mL of Anti-HA-tag Gel bound 1.7 mg of a HA-tagged protein (32 kDa) and eluted 1.4 mg of purified protein in our hands.

Storage

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

Material Preparation

Prepare the following reagents before affinity purification.

Lysis buffer : Suitable Lysis buffer varies with a kind of the HA-tagged protein
(see **Additional Information**).

10-50 mM Tris-HCl (pH 7.5)
100-300 mM NaCl
1-5 mM EDTA
1% NP-40 or Triton X-100
if necessary add Protease Inhibitor Cocktail
(e.g. SIGMA: code P8340, PIERCE: code 78415)


Washing buffer : PBS

Elution buffer : 2 mg/mL HA-tag peptide in PBS.

Regeneration buffer : 0.17 M Glycine-HCl, pH 2.3

Column storage buffer : 0.09% NaN₃/PBS

Procedure Summary

Protocols

- Notes:
1. The Gel is optimized under the native conditions only, and it is not recommended under denaturing conditions and also for purification of aggregated, unstable, and insoluble protein (e.g. inclusion bodies). Proteins solubilized with such as 6 M Guanidine-HCl or 8 M Urea can not be purified using this Gel (see **Additional Information**).
 2. Cellular debris and particulate matter must be removed prior to purification. The protein extract should be centrifuged (10,000-20,000 x g for 15 minutes) and filtered with a 0.45 µm filter to remove any remaining cells and particulates.
 3. Highly viscous samples containing chromosomal DNA or RNA should be sonicated or treated with nuclease to reduce viscosity.

A. Column preparation

1. Place the empty chromatography column (e.g., PIERCE, code 29920) on rack or stand.
2. Rinse the column with Washing buffer.
3. Resuspend the anti-HA-tag Gel by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
4. Transfer the desired volume to the column. Allow the column to drain. Do not allow the column to dry out.
5. Wash the column with 10 bed volumes of Washing buffer.

B. Column binding

1. Apply the lysate on the top of the column under gravity flow.

Note: The binding efficiency to anti-HA-tag Gel depends on the conditions such as a kind of the HA-tagged protein, sample flow rate, and temperature. If the coupling efficiency has been low, put a sample through the column several times, or incubate the lysate with the Gel in batch.

2. Collect the flow-through into clean collection tubes.
3. Wash the column with PBS until the OD280 is <0.01.

C. Elution of the HA-tagged protein

1. Prepare 5 bed volumes of Elution buffer.
2. Remove the bottom cap, and pour off the excess liquid.
3. To exchange the buffer in the column, add 1 bed volume of Elution buffer, and drain the 1 bed volume of the buffer.
4. Cap the bottom, and add 1 bed volumes of Elution buffer so that the Gel does not dry.

5. Allow the column to stand at 4°C for 5 minutes. (As the HA-tag peptide competes with HA-tagged protein on the Gel by this incubation, HA-tagged protein is eluted from the Gel.)
6. Remove the bottom cap, and collect the eluate into clean collection tubes (Fraction 1).
7. Add 1 bed volumes of Elution buffer, collect the eluate into clean collection tubes (Fraction 2).
8. Repeat the elution and collection step (Step 6 and 7) another 3 times. (Fraction 1-5 contains HA-tagged Protein.)

Note: Elution efficiency depends on each HA-tagged protein. If the efficiency has been low, please examine numbers of elution, and volume and peptide concentration of Elution buffer.

D. Regeneration and storage

1. Wash the column with 10 bed volumes of Regeneration buffer.
2. Immediately wash the column with 10 bed volumes of Column storage buffer.
3. Store at 2-8°C in 2 bed volumes of Column storage buffer.

Note: Poured columns containing the anti-HA-tag Gel may be used at least 10 times, depending on the usage conditions.


Related Products:

3320	HA-tagged Protein Purification Kit	20 purifications
3320A	HA-tagged Protein Purification Kit	2 purifications
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide	2 mg x 5
3325	DDDDK-tagged Protein Purification Kit	20 purifications
3325A	DDDDK-tagged Protein Purification Kit	2 purifications
3326	DDDDK-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 1, Peptide: 1 mg x 5
3327	DDDDK-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 5, Peptide: 1 mg x 25
3328	DDDDK-tagged Protein Purification Gel	5 mL x 1
3329	DDDDK-tagged Protein Purification Gel	5 mL x 5
3325-205	DDDDK-tag peptide	1 mg x 5
3310	His-tagged Protein Purification Kit	20 purifications
3310A	His-tagged Protein Purification Kit	2 purifications
3311	His-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 1, Peptide: 2 mg x 5
3312	His-tagged Protein Purification Gel with Elution Peptide	Gel: 1 mL x 5, Peptide: 2 mg x 25
3310-205	His-tag peptide	2 mg x 5

3305	c-Myc-tagged Protein Mild Purification Kit ver.2	20 purifications
3305A	c-Myc-tagged Protein Mild Purification Kit ver.2	2 purifications
3306	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3307	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 5	
3300-205	c-Myc-tag peptide	1 mg x 5
3315	V5-tagged Protein Purification Kit	20 purifications
3315A	V5-tagged Protein Purification Kit	2 purifications

Example of Purification Results

Purification of N-terminus HA-tagged β -galactosidase

SDS-PAGE (Coomassie Brilliant Blue Staining)

Human embryonic kidney cells (293T) were transfected with pcDNA-HA-tagged β -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (10 mL/100-mm dish x5) and purified according to the preceding protocol. Column bed volume was 0.25 mL. Elution was carried out with 1.5mL of 2 mg/mL HA-tag peptide. Each fraction was 0.25 mL.

Additional Information

Several reagents were examined whether or not they were suitable for use with the HA-tagged Protein Purification Gel. For example, RIPA buffer could be used for preparation of cell lysate. The results are listed below.

Chaotropic agents

Urea	1 M	Yes
Guanidine-HCl	1 M	No

Reducing agents

DTT	10 mM	Yes
2-Mercaptoethanol	10 mM	Yes

Surfactants

Nonionic	Tween-20	5%	Yes
	Triton X-100	5%	Yes
	NP40	5%	Yes
	Digitonin	1%	Yes
	n-Octyl-β-D-galcoside	1%	Yes
	CHAPS	1%	Yes
Zwitterionic	CHAPSO	1%	Yes
	SDS	0.05%	No
Anionic	Sodium Deoxycholate	0.5%	Yes

Others

NaCl	1 M	Yes
Glycerol	10%	Yes
EDTA	10 mM	Yes

The “Yes” indicates the reagents can be used in the Lysis buffer for this Gel up to the indicated concentration. The “No” indicates the reagents can not be used in the Lysis buffer for this Gel at the indicated concentration.

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
3320	HA-tagged Protein Purification Kit	20 purifications
3320A	HA-tagged Protein Purification Kit	2 purifications
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide	2 mg x 5
3325	DDDDK-tagged Protein Purification Kit	20 purifications
3325A	DDDDK-tagged Protein Purification Kit	2 purifications
3326	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3327	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 25	
3328	DDDDK-tagged Protein Purification Gel	5 mL x 1
3329	DDDDK-tagged Protein Purification Gel	5 mL x 5
3325-205	DDDDK-tag peptide	1 mg x 5
3310	His-tagged Protein Purification Kit	20 purifications
3310A	His-tagged Protein Purification Kit	2 purifications
3311	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 2 mg x 5	
3312	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 2 mg x 25	
3310-205	His-tag peptide	2 mg x 5

3305	c-Myc-tagged Protein Mild Purification Kit ver.2	20 purifications
3305A	c-Myc-tagged Protein Mild Purification Kit ver.2	2 purifications
3306	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3307	c-Myc-tagged Protein Mild Purification Gel	
	Gel: 1 mL x 5, Peptide: 1 mg x 5	
3300-205	c-Myc-tag peptide	1 mg x 5
3315	V5-tagged Protein Purification Kit	20 purifications
3315A	V5-tagged Protein Purification Kit	2 purifications

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