Product datasheet MON-RTU1003



## Polyclonal antibody against ACTH

Clone no. -

MONOSAN Ready To Use

Product name Polyclonal antibody against ACTH

**Host** Rabbit

**Applications** IHC-P

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

**Clonality** Polyclonal

Clone number -

Size 7 ml

**Concentration** n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON-RTU1003

**MONOSAN** 

Polyclonal antibody against ACTH

Clone no.

MONOSAN Ready To Use

## Additional info

ACTH or Adrenocorticotropic hormone is synthesized from pre-proopiomelanocortin (pre-POMC). ACTH is produced and secreted from corticotrophs in the anterior lobe (or adenohypophysis) of the pituitary gland. The anti-ACTH immunohistochemical reagent could be useful in the study of neoplastic and non-neoplastic pituitary diseases

References 1. Pizarro CB, et al. Braz J Med Biol Res. 2004; 37:235-43

2 Kageyama K, et al. Am J Med Sci. 2002; 324:326-30

3. Fan X, et al. J Histochem Cytochem. 2002; 50:1509-16

4. Japon MA, et al. J Clin Endocrinol Metab. 2002; 87:1879-84

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES