

Polyclonal antibody against CD3

Clone no. -

MONOSAN Ready To Use

---

Product name	Polyclonal antibody against CD3
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against CD3

Clone no. -

MONOSAN Ready To Use

---

**Additional info**

Anti-CD3 antibody has been considered the best all around T-cell marker. This antibody reacts with an antigen present in early thymocytes. The positive staining of this marker may represent a sign of early commitment to the T-cell lineage.

**References**

1. Beverley PC, et al. Eur J Immunol. 1981; 11:329-34
2. Clevers H, et al. Eur J Immunol. 1988; 18:705-10
3. Hedvat CV, et al. Hum Pathol. 2002; 33:968-74
4. Karube K, et al. Am J Surg Pathol. 2003; 27:1366-74
5. Dogan A, et al. Am J Surg Pathol. 2003; 27:903-11

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES