

Mouse anti-CD45RO, clone UCHL-1

Clone no. UCHL-1

MONOSAN Ready To Use

---

Product name	Mouse anti-CD45RO, clone UCHL-1
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a-k
Clonality	Monoclonal
Clone number	UCHL-1
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD45RO, clone UCHL-1

Clone no. UCHL-1

MONOSAN Ready To Use

---

**Additional info**

Anti-CD45RO labels an isoform of the CD45 antigen also known as leukocyte common antigen. Anti-CD45RO reacts with thymocytes, mature activated T-cells, and a subpopulation of resting T-cells while showing no reactivity with B-cells, making this antibody helpful in identifying T-cell neoplasms.

**References**

1. Hall PA, et al. J Clin Pathol. 1987; 40:151-6
2. Smith SH, et al. Immunology 1986; 58:63-70
3. Falini B, et al. Hum Pathol. 1990; 21:624-9
4. Koch AE, et al. J Clin Immunol. 1990; 10:192-9
5. Ritter JH, et al. J Cutan Pathol. 1994; 21:481-93

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES