

Polyclonal antibody against Gastrin

Clone no. -

MONOSAN Ready To Use

Product name	Polyclonal antibody against Gastrin
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against Gastrin

Clone no. -

MONOSAN Ready To Use

Additional info

Anti-Gastrin antibody gives positive staining of G-cells of human antral/pyloric mucosa and cells producing gastrin or a structural gastrin analogue as is seen in stomach; no staining of other cells or tissue types has been observed. This antibody may react with sulfated and non-sulfated forms of gastrin. The antibody cross-reacts with more than 50% of the present cholecystokinin octapeptide.

References

1. Kasacka W, et al. Folia Morphol. 2012; 71:39-44.
2. Hur K, et al. J Cancer Res Clin Oncol. 2006; 132:85-91
3. Waldum et al. Frontiers in Endocrinology. 2017; 8:1-7
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES