

## Polyclonal antibody against IgG

Clone no. -

MONOSAN Ready To Use

Product name Polyclonal antibody against IgG

**Host** Rabbit

**Applications** IHC-P

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

**Clonality** Polyclonal

Clone number -

Size 7 ml

**Concentration** n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON-RTU1131

**MONOSAN** 

Polyclonal antibody against IgG

Clone no. - MONOSAN Ready To Use

## Additional info

Anti-lgG reacts with immunoglobulin gamma (lgG) chains. This antibody is useful when differentiating and sub-classifying hematolymphoid neoplasms.

**References** 1. Arnold A, et al. New Eng J Med. 1983; 309:1593-1599

2 Leong AS, et al. Geenwich Medical Media Ltd. 1999; 217-219

3. Ando K, et al. Intern Med. 2000; 39:170-5

4. Taylor CR. Arch Path Lab Med. 1978; 102:113-121

5. Warnake R, et al. Masson Publishing USA. 1981; 203-221

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES