### Product datasheet

MON-RTU1141



MONOSAN Ready To Use

Polyclonal antibody against Lysozyme Clone no. -

Product name	Polyclonal antibody against Lysozyme
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	-
Clonality	Polyclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	_
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

### Product datasheet

MON-RTU1141

Polyclonal antibody against Lysozyme Clone no.

# **MONOSAN**<sup>®</sup>

MONOSAN Ready To Use

### Additional info

Anti-Lysozyme stains myeloid cells, histiocytes, granulocytes, macrophages, and monocytes in human tonsil, colon and skin. It is an important marker that may demonstrate the myeloid or monocytic nature of acute leukemia. The restrictive nature of anti-lysozyme antibody staining suggests that lysozyme may be synthesized predominantly in reactive histiocytes rather than in resting, unstimulated phagocytes. Anti-lysozyme may aid in the identification of histiocytic neoplasias, large lymphocytes and classifying lymphoproliferative disorders.

Rehg J, et al. Toxicol Pathol. 2012; 40: 345-74 Seifert RP, et al. Annals Diag Pathol. 2014; 18:253-60.

- 3. -
- 4. -
- 5. -

1. 2

#### FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com