

Mouse anti-Melanoma, clone HMB-45

Clone no. HMB-45

MONOSAN Ready To Use

Product name	Mouse anti-Melanoma, clone HMB-45
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1-k
Clonality	Monoclonal
Clone number	HMB-45
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Melanoma, clone HMB-45

Clone no. HMB-45

MONOSAN Ready To Use

Additional info

Anti-HMB45 is a useful melanoma immunohistochemical marker that reacts with antigens present on immature melanosomes. Anti-HMB45 is useful for identifying amelanotic melanoma from other neoplastic lesions with similar morphology.

References

1. Gown AM, et al. Am J Pathol. 1986; 123:195-203
2. Wick MR, et al. Arch Pathol Lab Med. 1988; 112:616-20
3. Abrahamsen HN, et al. Cancer. 2004; 100:1683-91
4. Vaggelli L, et al. Tumori. 2000; 86:346-8
5. Baisden BL, et al. Am J Surg Pathol. 2000; 24:1140-6

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES