

Polyclonal antibody against Myeloperoxidase

Clone no. -

MONOSAN Ready To Use

Product name	Polyclonal antibody against Myeloperoxidase
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against Myeloperoxidase

Clone no. -

MONOSAN Ready To Use

Additional info

Anti-Myeloperoxidase detects granulocytes and monocytes in blood and precursors of granulocytes in the bone marrow. This antibody can detect myeloid cell populations of the bone marrow as well as in other sites.

References

1. Pinkus GS, et al. Mod Pathol. 1991; 4:733-41
2. Markoc F, et al. Tumori. 2010; 96:149-53
3. Alexiev BA, et al. Diagn Pathol. 2007; 31;2:42
4. Saravanan L, et al. Int J Lab Hematol. 2010; 32:132-6
5. Manaloor EJ, et al. Am J Clin Pathol. 2000; 113:814-22

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES