

Mouse anti-PMS2, clone MRQ-28

Clone no. MRQ-28

MONOSAN Ready To Use

---

Product name	Mouse anti-PMS2, clone MRQ-28
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	MRQ-28
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-PMS2, clone MRQ-28

Clone no. MRQ-28

MONOSAN Ready To Use

---

**Additional info**

Analyte Specific Reagent: Analytical and performance characteristics are not determined.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES