

Polyclonal antibody against Prolactin

Clone no. -

MONOSAN Ready To Use

Product name	Polyclonal antibody against Prolactin
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against Prolactin

Clone no. -

MONOSAN Ready To Use

Additional info

Prolactin (PRL) is a single-chain polypeptide of 226 amino acids and plays a role in multiple processes including cell growth, reproduction, and immune function. Anti-Prolactin reacts with prolactin-producing cells and is a useful marker in classification of pituitary tumors and the study of pituitary disease.

References

1. Asa SL, et al. Arch Pathol Lab Med. 1982; 106:360-3
2. Duello TM, et al. Am J Anat. 1980; 158:463-9
3. Minniti G, et al. Surg Neurol. 2002; 57:99-103
4. Popadic A, et al. Surg Neurol. 1999; 51:47-54
5. Nevalainen MT, et al. J Clin Invest. 1997; 99:618-27

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES