

Mouse anti-Spectrin, clone RBC2/3D5

Clone no. RBC2/3D5

MONOSAN Ready To Use

| | |
|---------------------------|---|
| Product name | Mouse anti-Spectrin, clone RBC2/3D5 |
| Host | Mouse |
| Applications | IHC-P |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | Unknown or proprietary to MONOSAN and/or its suppliers |
| Isotype | IgG2b-k |
| Clonality | Monoclonal |
| Clone number | RBC2/3D5 |
| Size | 7 ml |
| Concentration | n/a |
| Format | - |
| Storage buffer | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Spectrin, clone RBC2/3D5

Clone no. RBC2/3D5

MONOSAN Ready To Use

Additional info

Spectrin is a cytoskeletal protein which is found in muscles, red blood cells and red cell precursors. Anti-Spectrin antibody is useful in the identification of blood dyscrasias and muscle disorders.

References

1. Sadahira, Y et al. J Clin Pathol. 1999 Dec; 52(12): 919-21
2. Nehls, V et al. Am J Pathol. 1993 May; 142(5): 1565-73
3. Muller M et al. J Vet Med A Physiol Pathol Clin Med. 2001 Feb;48(1):51-7
4. Terada N et al. J Anat. 1997 Apr;190(Pt 3):397-404
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES