### Product datasheet

MON-RTU1195



Polyclonal antibody against Synaptophysin Clone no. -

		_	
MONOSAN	Ready	TO	I JSP
111011000/ 111	ricaay	. 0	000

Product name	Polyclonal antibody against Synaptophysin
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	-
Clonality	Monoclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

### Product datasheet

MON-RTU1195

Polyclonal antibody against Synaptophysin Clone no. -

# MONOSAN

MONOSAN Ready To Use

### Additional info

Anti-Synaptophysin reacts with neuroendocrine cells of human adrenal medulla, carotid body, skin, pituitary, thyroid, lung, pancreas and gastrointestinal mucosa. Positive staining is seen in neurons of the brain, spinal cord, retina, and Paneth's cells in the gastrointestinal tract and gastric parietal cells. This antibody identifies normal neuroendocrine cells and neuroendocrine neoplasms. Diffuse, finely granular cytoplasmic staining is observed, which probably correlates with the distribution of the antigen within neurosecretory vesicles. The expression of synaptophysin is independent of the presence of NSE or other neuroendocrine markers. Anti-Synaptophysin is an independent broadrange marker of neural and neuroendocrine differentiation.

#### References

1.

Navone F, et al. J Cell Biol. 1986; 103:2511-27

- 2 Wiedenmann B, et al. Cell. 1985; 41:1017-28
- 3. Kayser K, et al. Pathol Res Pract. 1988; 183:412-7
- 4. Conner MG, et al.Ann Diagn Pathol. 2002; 6:345-8
- 5. Lyda MH, et al. Hum Pathol. 2000; 31:980-7

### FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com