Product datasheet MON-RTU1200



Polyclonal antibody against Toxoplasma gondii

Clone no.

MONOSAN Ready To Use

Product name Polyclonal antibody against Toxoplasma gondii

Host Rabbit

Applications IHC-P

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

Clonality Monoclonal

Clone number -

Size 7 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON-RTU1200



Polyclonal antibody against Toxoplasma gondii

Clone no.

MONOSAN Ready To Use

Additional info

Toxoplasma gondii is a spindle-to-oval-shaped protozoan which presents as an infection in humans of various sorts. The cyst (30 um) and trophozoite (7 um) stages can be identified in humans is such cases. This intracellular parasite is transmitted via raw/undercooked meat, contaminated soil, or by direct contact with an infected host. Infection in humans is usually associated with a variable degree of immunosuppression such as in pregnancy or immunosuppression due to various drugs.

References 1. -

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES