

Mouse anti-MART-1, clone A103

Clone no. A103

MONOSAN Ready To Use

Product name	Mouse anti-MART-1, clone A103
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	A103
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MART-1, clone A103

Clone no. A103

MONOSAN Ready To Use

Additional info

MART-1 (also known as Melan A) is a melanocyte differentiation antigen. MART-1 is a transmembrane protein present in melanocytes of normal skin, retina, nevi, and most melanomas. MART-1 is a very useful marker for identifying metastatic melanomas.

References

1. Kageshita T et al. J Immunother 1997 Nov;20(6):460-5
2. Yaziji H, et al. In J Surg Pathol. 2003 Jan;11(1):11-5
3. Mocellin S et al. J Immunother. 2001 Nov-Dec;24(6):447-58
4. Perez RP et al. Hum Pathol. 2000 Nov;31(11):1381-8
5. Hoang MP et al. J Cutan Pathol. 2001 Sep;28(8):400-6

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES