

Mouse anti-ksp-Cadherin, clone MRQ-33

Clone no. MRQ-33

MONOSAN Ready To Use

---

Product name	Mouse anti-ksp-Cadherin, clone MRQ-33
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	MRQ-33
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-ksp-Cadherin, clone MRQ-33

Clone no. MRQ-33

MONOSAN Ready To Use

**Additional info**

Kidney-specific cadherin (Ksp-cadherin) is a member of the cadherin family of cell adhesion molecules that is found exclusively in the basolateral membrane of renal tubular epithelial cells of the distal tubules and collecting duct. Ksp-cadherin may be useful in distinguishing between renal neoplasms of distal nephron origin from those of proximal tubule origin.

**References**

1. Mazal PR, et al. Hum Pathol. 2005; 36:22-8
2. Shen SS, et al. Mod Pathol. 2005; 18:933-40
3. Thedieck C, et al. Br J Cancer. 2005; 92:2010
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES