

Rabbit anti-CD3, clone MRQ-39

Clone no. MRQ-39

MONOSAN Ready To Use

---

Product name	Rabbit anti-CD3, clone MRQ-39
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG
Clonality	Monoclonal
Clone number	MRQ-39
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rabbit anti-CD3, clone MRQ-39

Clone no. MRQ-39

MONOSAN Ready To Use

---

**Additional info**

CD3's immunohistochemical detection locates the cytoplasmic component of CD3 protein. Anti-CD3 is considered to be a pan-T-cell marker and reacts with an antigen present at the surface and in the cytoplasm of T lymphocytes. Anti-CD3 is widely used for the identification of immature and mature T-cell malignancies.

**References**

1. Beverley PC, et al. Eur J Immunol.; 11:329-34 (1981)
2. Hedvat CV, et al. Hum Pathol.; 33:968-74 (2002)
3. Karube K, et al. Am J Surg Pathol.; 27:1366-74 (2003)
4. Dogan A, et al. Am J Surg Pathol.; 27:903-11 (2003)
5. Axdorph U, et al. APMIS.; 110:379-90 (2002).

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**