

Mouse anti-Melanoma Associated Antigen, clone KBA.62

Clone no. KBA.62

MONOSAN Ready To Use

Product name	Mouse anti-Melanoma Associated Antigen, clone KBA.62
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	KBA.62
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Melanoma Associated Antigen, clone KBA.62

Clone no. KBA.62

MONOSAN Ready To Use

Additional info

The antibody is a monoclonal, anti-melanoma antibody that reacts with an antigen that has yet to be identified.¹ Notably used as a melanoma marker, KBA.62 also detects smooth muscle, basal cells of the epidermis and hair shaft epithelia of the skin.

References

1. Pages C, et al.. Hum Pathol. 2008; 39:1136-42
2. Aung P, et al. Am J Surg Pathol. 2012; 36:265-72
3. E Cohen-Knafo et al. J Clin Pathol. 1995; 48:826-831
4. Kaufmann O, et al.. Mod Pathol, 1998 Aug; 11(8):740-6
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES