

## Mouse anti-MiTF, clone C5/D5

Clone no. C5/D5

MONOSAN Ready To Use

---

Product name	Mouse anti-MiTF, clone C5/D5
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	C5/D5
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MiTF, clone C5/D5

Clone no. C5/D5

MONOSAN Ready To Use

---

**Additional info**

Mi is a transcription factor implicated in pigmentation, bone development and in mast cells. Various forms of Mi exist ranging from 50-70 kD in size. This antibody targets the 52-56 kD range. This antibody has been useful in identifying malignant melanoma.

**References**

1. Liegl B, et al. Am J Surg Pathol. 2008; 32:608-14
2. Righi A, et al. Int J Surg Pathol. 2008; 16:16-20
3. Weinreb I, et al.. Virchows Arch. 2007; 450:463-70
4. Ohsie SJ, et al. J Cutan Pathol. 2008; 35:433-44
5. Hornick JL, et al. Am J Surg Pathol. 2008; 32:493-501

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES