

Mouse anti-Cytokeratin, clone OSCAR

Clone no. OSCAR

MONOSAN Ready To Use

Product name	Mouse anti-Cytokeratin, clone OSCAR
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a
Clonality	Monoclonal
Clone number	OSCAR
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Cytokeratin, clone OSCAR

Clone no. OSCAR

MONOSAN Ready To Use

Additional info

Anti-Cytokeratin (OSCAR) is well suited to distinguish epithelial carcinoma from non-epithelial malignancies and is used to aid epithelial tumor classification. This antibody has been used to characterize the source of various neoplasms and to study the distribution of keratin containing cells in epithelia during normal development and during the development of epithelial neoplasms. This antibody stains cytokeratins present in normal and abnormal human tissues and has shown high sensitivity in recognizing epithelial cells, and carcinomas.

References

1. Gown, AM, et al. Am J Clin Pathol 1985;84:413
2. Battifora, H. Am J Surg Pathol 1988;12:24
3. Lewis JE et al. Hum Pathol. 1997 Jun;28(6):664-73
4. Mueller JD et al. Cancer. 2000 Nov 1;89(9):1874-82
5. Sato F et al. Br J Surg. 2001 Mar;88(3):426-32

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES