Product datasheet MON-RTU1284



Mouse anti-Hepatitis B Virus Surface Antigen, clone A10F1

Clone no. A10F1 MONOSAN Ready To Use

Product name Mouse anti-Hepatitis B Virus Surface Antigen, clone A10F1

**Host** Mouse

**Applications** IHC-P

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

**Clonality** Monoclonal

Clone number A10F1

Size 7 ml

**Concentration** n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON-RTU1284



Mouse anti-Hepatitis B Virus Surface Antigen, clone A10F1

Clone no. A10F1

MONOSAN Ready To Use

Additional info

References 1. -

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES