Product datasheet MON-RTU1290



Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Clone no. DBA.44 MONOSAN Ready To Use

Product name Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Host Mouse

Applications IHC-P (1:15-1:30)

Species reactivity human

Conjugate -

Immunogen DEAU cell line established from a large cell lymphoma

Isotype IgMk

Clonality Monoclonal

Clone number DBA.44

Size 6 ml

Concentration lot specific

Format -

Storage buffer Tissue culture supernatant with sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON-RTU1290



Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Clone no. DBA.44 MONOSAN Ready To Use

Additional info

This antibody reacts with a subset of B-lymphocytes localized in the follicular mantle zone. It reacts with 97% of hairy cell leukemia cases. This antibody shows strong positive staining of about 35% of cases of high grade B cell lymphomas. Positive control Tonsil.

References 1. Al Saati et al. 1989;74:2476

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES