

Mouse anti-CD8, clone FK18 (Monoclonal)

Clone no. FK18

MONOSAN

---

Product name	Mouse anti-CD8, clone FK18 (Monoclonal)
Host	Mouse
Applications	IF (1:10), preferably in PBS plus 1% BSA and 0.02% sodium azide
Species reactivity	human, rhesus
Conjugate	-
Immunogen	Human thymocytes
Isotype	IgG3
Clonality	Monoclonal
Clone number	FK18
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Serum free medium containing 1% FCS
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-CD8, clone FK18 (Monoclonal)

Clone no. FK18

MONOSAN

**Additional info**

The antibody reacts with CD8, the 33 kD T8 antigen present on the human cytotoxic/suppressor T-cell subset. It strongly blocks the cytolytic activity of T8-positive CTLs. It also reacts with a subpopulation of Rhesus monkey lymphocytes which is identical to the population of human T8 positive cells.

**References**

1. Jansen, J., et al., 1984, Seminars in Hematology 21, 164
2. Spits, H., et al., 1983, Hybrydoma 2, 423
3. Roelen, D.L., et al., 1992, Transduction 54, 899
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES