

Mouse anti-CD57, clone 6-13-19-1 (Monoclonal)

Clone no. 6-13-19-1

MONOSAN

Product name	Mouse anti-CD57, clone 6-13-19-1 (Monoclonal)
Host	Mouse
Applications	IF (1:5), IHC (1:5)
Species reactivity	human
Conjugate	-
Immunogen	JURKAT T-cell line
Isotype	IgM
Clonality	Monoclonal
Clone number	6-13-19-1
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD57, clone 6-13-19-1 (Monoclonal)

Clone no. 6-13-19-1

MONOSAN

Additional info

The antibody reacts with a 110 kD antigen present on the peripheral blood lymphocyte subset (primarily large granular lymphocytes) identified by the HNK-1 monoclonal antibody. It recognizes part of the Natural Killer cell population. CD57 is furthermore present on cells in neuroectodermal tissue

References

1. Abo, T., and Balch, C.M., 1981, The Journal of Immunology 127, 1024
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES