### Product datasheet

MON1027



Mouse anti-CD15, clone BRA4F1 (Monoclonal) Clone no. BRA4F1

Product name	Mouse anti-CD15, clone BRA4F1 (Monoclonal)
Host	Mouse
Applications	ELISA, FC, IHC-fr, IHC-P, IF
Species reactivity	human
Conjugate	-
Immunogen	erythroid-myeloid leukemia cell line
lsotype	IgM
Clonality	Monoclonal
Clone number	BRA4F1
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C
-0	

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Product datasheet

MON1027

# MONOSAN

Mouse anti-CD15, clone BRA4F1 (Monoclonal)Clone no.BRA4F1

MONOSAN

### Additional info

Bra4F1 reacts with CD15 (220 kDa). CD15 is present on >95% of granulocytes including neutrophils and eosinophils and to a lesser degree on monocytes. CD15 is further expressed in Reed-Sternberg cells in classic Hodgkin's disease. CD15 is occasionally expressed in large cell lymphomas of both B and T phenotypes. It is also expressed on a wide variety of other tumor cells including myeloid leukemia, breast, colorectal, and lung cancer cells. Bra4F1 was clustered at the IVth International Workshop on Leucocyte Differentiation Antigens.

References	1.	Chorvath B. et al., Neoplasma 36(6), 633-642 (1989)
	2	Leukocyte typing IV, Oxford Univ. Press, 868, 870 and 877 (1989)
	3.	-
	4.	-
	5.	-

### FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021