

Product datasheet MON10302

MONOSAN[®]

Mouse anti-Human Aggrecan, clone 70 (Monoclonal)

Clone no. 70

MONOSAN

Product name Mouse anti-Human Aggrecan, clone 70 (Monoclonal)

Host Mouse

Applications IHC-fr, ELISA, IP, WB

Species reactivity human, bovine

Conjugate -

Immunogen Aggrecan

Isotype IgG1

Clonality Monoclonal

Clone number 70

Size 0.2 mg

Concentration IgG 1 mg/ml

Format Purified

Storage buffer PBS with <0.1% sodium azide

Storage until expiry date aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021

Revision date
16-11-2022

Mouse anti-Human Aggrecan, clone 70 (Monoclonal)

Clone no.

70

MONOSAN

Additional info

Mouse anti Human aggrecan antibody, clone 7E1 recognizes human aggrecan, a proteoglycan and member of the aggrecan/versican proteoglycan family, which forms a major component of the extracellular matrix (ECM) of both cartilage and the central nervous system (CNS).

References

1. Virgintino D et al. J Cell Mol Med 2009; 13: 3151-73
2. Aigner T et al. Cancer 2002; 94: 2273-81
3. Wei A et al. J Orthop Res 2016; 34: 279-89
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES