

Mouse anti-Human Alkaline Phosphatase, clone BGN/03/661 (Monoclonal)

Clone no. BGN/03/661

MONOSAN

Product name	Mouse anti-Human Alkaline Phosphatase, clone BGN/03/661 (Monoclonal)
Host	Mouse
Applications	ELISA, IHC-P, WB
Species reactivity	human
Conjugate	-
Immunogen	Alkaline phosphatase from human bone
Isotype	IgG1
Clonality	Monoclonal
Clone number	BGN/03/661
Size	0.2 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Alkaline Phosphatase, clone BGN/03/661 (Monoclonal)

Clone no. BGN/03/661

MONOSAN

Additional info

Mouse anti Human Alkaline Phosphatase antibody, clone BGN/03/661 recognizes human alkaline phosphatase, a hydrolase enzyme most active at alkaline pH. Mouse anti Human Alkaline Phosphatase antibody, clone BGN/03/661 detects the tissue specific isotypes, BAP, PLAP and IAP (bone, placental and intestinal alkaline phosphatase), but may react with other alkaline phosphatases.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES