

Mouse anti-Human Beta Endorphin, clone B31.15 (Monoclonal)

Clone no. B31.15

MONOSAN

Product name	Mouse anti-Human Beta Endorphin, clone B31.15 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IF, IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Human beta endorphin (full length native protein)
Isotype	IgG1
Clonality	Monoclonal
Clone number	B31.15
Size	0.2 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Beta Endorphin, clone B31.15 (Monoclonal)

Clone no. B31.15

MONOSAN

Additional info

Mouse anti Human beta endorphin antibody, clone B31.15 reacts with beta endorphin, a 31 amino acid, opioid neurotransmitter found in the neurons of the central and peripheral nervous system

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES