

Mouse anti-Human CD79a, clone HM57 (Monoclonal)

Clone no. HM57

MONOSAN

Product name	Mouse anti-Human CD79a, clone HM57 (Monoclonal)
Host	Mouse
Applications	IHC-fr, FC, IHC-P, WB
Species reactivity	human, mouse, rabbit, rat, pig, bovine, guinea pig, monkey
Conjugate	-
Immunogen	Synthetic peptide corresponding to 202-216 amino acid sequence of human mb-1
Isotype	IgG1
Clonality	Monoclonal
Clone number	HM57
Size	0.2 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human CD79a, clone HM57 (Monoclonal)

Clone no. HM57

MONOSAN

Additional info

Mouse anti Human CD79a antibody, clone HM57 recognizes an epitope within the cytoplasmic domain of CD79a. CD79a, also known as mb-1, is a 45 kDa protein that is expressed by B lymphocytes during differentiation from early pre-B cell stage through to plasma cells. The CD79a molecule associates with CD79b (B29) to form a heterodimer that is non-covalently linked to surface immunoglobulin, forming the B-cell receptor (BCR) complex. The CD79a/CD79b heterodimers are also necessary for intracellular signaling following antigen-binding to surface immunoglobulin.

References

1. Mason DY et al. J Immunol 1991; 147: 2474-82
2. Spaas JH Vet J 2013; 195: 107-13
3. De Schauwer C et al. Stem Cell Res Ther 2014; 5: 6
4. Hillmann A et al. Cell transplant 2016; 25: 109-24
5. Uitterdijk A et al. PLoS One 2017; 12: e0178779

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES