

Mouse anti-Human Insulin/Proinsulin, clone D3E7 (5B6/6) (Monoclonal)

Clone no. D3E7 (5B6/6)

MONOSAN

---

Product name	Mouse anti-Human Insulin/Proinsulin, clone D3E7 (5B6/6) (Monoclonal)
Host	Mouse
Applications	IHC-fr, ELISA, IHC-P
Species reactivity	human, mouse, rat, bovine, pig
Conjugate	-
Immunogen	Recombinant human insulin
Isotype	IgG1
Clonality	Monoclonal
Clone number	D3E7 (5B6/6)
Size	0.2 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Insulin/Proinsulin, clone D3E7 (5B6/6) (Monoclonal)

Clone no. D3E7 (5B6/6)

MONOSAN

**Additional info**

Mouse anti Human Insulin/Proinsulin antibody, clone D3E7 (5B6/6) recognizes both insulin and proinsulin, but does not react with free C-peptide.

**References**

1. Qin J et al. Life Sci 2009; 84: 810-6
2. Jörns A et al. Diabetes 2015; 64: 2880-91
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**