

Mouse anti-Human Lysozyme, clone BGN/0696/5B1 (Monoclonal)

Clone no. BGN/0696/5B1

MONOSAN

Product name	Mouse anti-Human Lysozyme, clone BGN/0696/5B1 (Monoclonal)
Host	Mouse
Applications	ELISA, FC, IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Native, from urine
Isotype	IgG2a
Clonality	Monoclonal
Clone number	BGN/0696/5B1
Size	0.2 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Lysozyme, clone BGN/0696/5B1 (Monoclonal)

Clone no. BGN/0696/5B1

MONOSAN

Additional info

Mouse anti Human Lysozyme antibody, clone BGN/0696/5B1 recognizes lysozyme, an enzyme that is present in mucosal secretions such as saliva and tears and is involved with non-specific organism defence. It functions by hydrolyzing 1,2-beta linkages between N-acetylmuramic acid and N-acetyl-D-glucosamine residues such as those found in the walls of various Gram-positive bacteria.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES