

Mouse anti-Human Troponin I (CARDIAC), clone 10F4 (110) (Monoclonal)

Clone no. 10F4 (110)

MONOSAN

Product name	Mouse anti-Human Troponin I (CARDIAC), clone 10F4 (110) (Monoclonal)
Host	Mouse
Applications	ELISA, IHC-P, WB
Species reactivity	human, rat, goat, cat, rabbit, bovine, dog, pig
Conjugate	-
Immunogen	Human Cardiac troponin I.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	10F4 (110)
Size	0.25 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Troponin I (CARDIAC), clone 10F4 (110) (Monoclonal)

Clone no. 10F4 (110)

MONOSAN

Additional info

Mouse anti Human Troponin I (Cardiac) antibody, clone 10F4 (110) recognizes aa 34-37 of Troponin I from human heart. There is no cross-reactivity with skeletal muscle troponin I.

References

1. Marian AJ et al. Circ Res 1997; 81: 76-85
2. Zhang J et al. Toxicol Pathol 2008; 36: 277-8
3. Filatov VL et al. Biochem Mol Biol Int 1998; 45: 1179-87
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES