

Product datasheet MON10400

**MONOSAN**<sup>®</sup>

Rabbit anti-Lysozyme

Clone no. -

MONOSAN

---

Product name	Rabbit anti-Lysozyme
Host	Rabbit
Applications	IHC-P (1:50-1:100)
Species reactivity	human
Conjugate	-
Immunogen	Lysozyme isolated from urine of patients with monocytic leukemia
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	Purified
Storage buffer	Immunoglobulin fraction of rabbit antiserum and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

**Additional info**

This antibody reacts with lysozyme. It stains granulocytes, monocytes and macrophages in human tonsil, skin and colon. This antibody does not crossreact with any other cell types.

**References**

1. Morsky et al. Clin. Chim. Acta 1988; 178: 327
2. Krugliak et al. Am J Hematol 1986; 21: 99
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES