

Mouse anti-HLA-DR+DP, clone HL-40 (Monoclonal)

Clone no. HL-40

MONOSAN

Product name	Mouse anti-HLA-DR+DP, clone HL-40 (Monoclonal)
Host	Mouse
Applications	FC , IP, IHC-P
Species reactivity	Human, Pig
Conjugate	-
Immunogen	Raji Burkitt's lymphoma cell line
Isotype	IgG2a
Clonality	Monoclonal
Clone number	HL-40
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON1083

MONOSAN[®]

Mouse anti-HLA-DR+DP, clone HL-40 (Monoclonal)

Clone no. HL-40

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES
