

Mouse anti-CD15, clone V1MC6, R-PE (Monoclonal)

Clone no. V1MC6

MONOSAN

Product name	Mouse anti-CD15, clone V1MC6, R-PE (Monoclonal)
Host	Mouse
Applications	FC
Species reactivity	Human
Conjugate	R-PE
Immunogen	Human CD15
Isotype	IgM
Clonality	Monoclonal
Clone number	V1MC6
Size	0.5 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 4mg/ml BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD15, clone V1MC6, R-PE (Monoclonal)

Clone no. V1MC6

MONOSAN

Additional info

CD15 (Lewis X, Le(x); stage specific embryonic antigen-1, SSEA-1) is a trisaccharide determinant (3-fucosyl-N-acetyllactosamine) expressed on several glycolipids, glycoproteins and proteoglycans of various cell types, e.g. granulocytes, mast cells, monocytes, macrophages, cells of gastric mucosa, nervous system or various tumour cells.

References

1. Okoye-Okafor UC et al. Nat Chem Biol 2015
2. Yang Z et al. Science 1998
3. Barcena A et al. Transfusion 2011
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES