

Mouse anti-PECAM-1/CD31, clone BV8 (Monoclonal)

Clone no. BV8

MONOSAN

Product name	Mouse anti-PECAM-1/CD31, clone BV8 (Monoclonal)
Host	Mouse
Applications	FC,ELISA,IF,IP
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	BV8
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-PECAM-1/CD31, clone BV8 (Monoclonal)

Clone no. BV8

MONOSAN

Additional info

This monoclonal antibody binds to human PECAM-1 (platelet endothelial cell adhesion molecule-1; CD31) a specific component of endothelial cell junctions. PECAM-1 is also expressed in platelets and leukocytes.

References

1. Rival; Y et al. J Immunol 1996; 157: 1233
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES