

Mouse anti-HLA-Class I, clone W6/32, FITC (Monoclonal)

Clone no. W6/32

MONOSAN

Product name	Mouse anti-HLA-Class I, clone W6/32, FITC (Monoclonal)
Host	Mouse
Applications	FC
Species reactivity	Cow, Non-human primates, Human, Cat
Conjugate	FITC
Immunogen	Membrane of human tonsil cells
Isotype	IgG2a
Clonality	Monoclonal
Clone number	W6/32
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-HLA-Class I, clone W6/32, FITC (Monoclonal)

Clone no. W6/32

MONOSAN

Additional info

HLA-class I major histocompatibility (MHC) antigens are intrinsic membrane glycoproteins expressed on nucleated cells and noncovalently associated with an invariant beta2 microglobulin. They carry foreign determinants important for immune recognition by cytotoxic T cells, thus important for anti-viral and anti-tumour defence. Human HLA-class I antigens are represented by HLA-A, HLA-B and HLA-C molecules.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES