

Mouse anti-HLA-Class I, clone W6/32, PE (Monoclonal)

Clone no. W6/32

MONOSAN

---

Product name	Mouse anti-HLA-Class I, clone W6/32, PE (Monoclonal)
Host	Mouse
Applications	FC
Species reactivity	Cow, Non-human primates, Human, Cat
Conjugate	PE
Immunogen	Membrane of human tonsil cells
Isotype	IgG2a
Clonality	Monoclonal
Clone number	W6/32
Size	0.1 mg
Concentration	0.1 mg/ml
Format	-
Storage buffer	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-HLA-Class I, clone W6/32, PE (Monoclonal)

Clone no. W6/32

MONOSAN

**Additional info**

HLA-class I major histocompatibility (MHC) antigens are intrinsic membrane glycoproteins expressed on nucleated cells and noncovalently associated with an invariant beta2 microglobulin. They carry foreign determinants important for immune recognition by cytotoxic T cells, thus important for anti-viral and anti-tumour defence. Human HLA-class I antigens are represented by HLA-A, HLA-B and HLA-C molecules.

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**