

Mouse anti-beta-HCG, clone 2092 (Monoclonal)

Clone no. 2092

MONOSAN

Product name	Mouse anti-beta-HCG, clone 2092 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:10), IHC-fr (1:10), WB, ELISA
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	2092
Size	1 ml
Concentration	5 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON2001-1

MONOSAN[®]

Mouse anti-beta-HCG, clone 2092 (Monoclonal)

Clone no. 2092

MONOSAN

Additional info

The antibody does stain HCG containing human cells

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES