

Mouse anti-Campylobacter jejuni, clone 57-24 (Monoclonal)

Clone no. 57-24

MONOSAN

---

Product name	Mouse anti-Campylobacter jejuni, clone 57-24 (Monoclonal)
Host	Mouse
Applications	ELISA
Species reactivity	bacteria
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG3
Clonality	Monoclonal
Clone number	57-24
Size	1 ml
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide and 0.1% BSA
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Campylobacter jejuni, clone 57-24 (Monoclonal)

Clone no. 57-24

MONOSAN

**Additional info**

The antibody reacts with a cell wall component of *C. jejuni* and *C. coli*. The Hippurinate reaction was considered as a marker for *C. jejuni*. The antibody does not react with *C. fetus* or *C. lari*.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES