

Mouse anti-SLPI, clone 31, Biotin (Monoclonal)

Clone no.      31

MONOSAN

---

Product name	Mouse anti-SLPI, clone 31, Biotin (Monoclonal)
Host	Mouse
Applications	-
Species reactivity	-
Conjugate	Biotin
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	31
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet                      MON2005B



Mouse anti-SLPI, clone 31, Biotin (Monoclonal)

Clone no.                      31

MONOSAN

Additional info                      -

References                      1.                      -  
   2.                      -  
   3.                      -  
   4.                      -  
   5.                      -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES