Product datasheet MON2008-1



Mouse anti-S. epidermis, clone 17-5 (Monoclonal)

Clone no. 17-5 MONOSAN

Product name Mouse anti-S. epidermis, clone 17-5 (Monoclonal)

Host Mouse

Applications ELISA (1:10)

Species reactivity bacteria

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

lsotype lgG3

Clonality Monoclonal

Clone number 17-5

Size 1 ml

Concentration 200 ug/ml

Format -

Storage buffer RPMI with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON2008-1



Mouse anti-S. epidermis, clone 17-5 (Monoclonal)

Clone no. 17-5 MONOSAN

Additional info

The antibody reacts with several Staphylococcus strains using ELISA. It does not react with several members of the Enterobacteriaceae and not with S. Aureus. It is not tested for IHC-P. The epitope is not determined, but it might be a component of the cell wall. The antibody recognizes whole bacteria.

References 1. -

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES