

Mouse anti-Pseudomonas, clone 69-5

Clone no. 69-5

MONOSAN

Product name	Mouse anti-Pseudomonas, clone 69-5
Host	Mouse
Applications	ELISA, Agglutination test
Species reactivity	bacteria
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a
Clonality	Monoclonal
Clone number	69-5
Size	1 ml
Concentration	25 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide and 0.1% BSA
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Pseudomonas, clone 69-5

Clone no. 69-5

MONOSAN

Additional info

The antibody reacts with flagella of approximately 30% of *Pseudomonas aeruginosa* isolates. The antibody can be used for agglutination test.

References

1. Widjoatmodjo, M.N., et al., 1991, Eur. J. Clin. Microbio. Infect. Dis. 10, 935-38
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES