

Mouse anti-Proteus spp, clone 31-17 (Monoclonal)

Clone no. 31-17

MONOSAN

---

Product name	Mouse anti-Proteus spp, clone 31-17 (Monoclonal)
Host	Mouse
Applications	ELISA
Species reactivity	bacteria
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a
Clonality	Monoclonal
Clone number	31-17
Size	1 ml
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide and 0.1% BSA
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Proteus spp, clone 31-17 (Monoclonal)

Clone no. 31-17

MONOSAN

**Additional info**

The antibody reacts with bacteria in suspension and is therefore useful for detection of Proteus spp. in bacterial suspensions. Various systems can be used with this antibody to detect the presence of Proteus spp.. The antibody reacts with approximately 50% of the Proteus spp. isolates tested.

**References**

1. Widjoatmodjo, M.N., et al., 1991, Eur. J. Clin. Microbio. Infect. Dis. 10, 935-38
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES