

Mouse anti-H-FABP, clone 67D3, Biotin (Monoclonal)

Clone no. 67D3

MONOSAN

---

Product name	Mouse anti-H-FABP, clone 67D3, Biotin (Monoclonal)
Host	Mouse
Applications	-
Species reactivity	-
Conjugate	Biotin
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	67D3
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON2020B

**MONOSAN**<sup>®</sup>

Mouse anti-H-FABP, clone 67D3, Biotin (Monoclonal)

Clone no. 67D3

MONOSAN

---

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES