

Mouse anti-Rat CD4-antigen, clone 15-8A2

Clone no. 15-8A2

MONOSAN

---

Product name	Mouse anti-Rat CD4-antigen, clone 15-8A2
Host	Mouse
Applications	IHC-fr,FC
Species reactivity	rat
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	15-8A2
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Rat CD4-antigen, clone 15-8A2

Clone no. 15-8A2

MONOSAN

**Additional info**

The monoclonal antibody recognizes rat CD4 antigen (MW 52 kDa) and is reactive with all common inbred strains. 15-8A2 is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from Balb/c mice immunized with WAG/Rij spleen cells.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES