

Mouse anti-Rat CD8, clone 15-11C5 (Monoclonal)

Clone no. 15-11C5

MONOSAN

Product name	Mouse anti-Rat CD8, clone 15-11C5 (Monoclonal)
Host	Mouse
Applications	IHC-fr,FC
Species reactivity	rat
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	15-11C5
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Rat CD8, clone 15-11C5 (Monoclonal)

Clone no. 15-11C5

MONOSAN

Additional info

The monoclonal antibody recognizes rat CD8 antigen (MW 34 kDa and 39 kDa) and is reactive with all common inbred strains. 15-11C5 is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from Balb/c mice immunized with WAG/Rij spleen cells.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES