

Mouse anti-Human blood group antigen A1B, clone HE-24 (Monoclonal)

Clone no. HE-24

MONOSAN

Product name	Mouse anti-Human blood group antigen A1B, clone HE-24 (Monoclonal)
Host	Mouse
Applications	FUNC
Species reactivity	Human
Conjugate	-
Immunogen	Glycoprotein fraction isolated from saliva of secretors with blood group A.
Isotype	IgM
Clonality	Monoclonal
Clone number	HE-24
Size	1.0 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with <0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human blood group antigen A1B, clone HE-24 (Monoclonal)

Clone no. HE-24

MONOSAN

Additional info

The antibody HE-24 distinguishes blood group A1B from A2B. The specificity of the antibody HE-24 was verified on >1000 samples of blood. Agglutination test A1, A1B.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES