

Mouse anti-Human blood group antigen B, clone HEB-29 (Monoclonal)

Clone no. HEB-29

MONOSAN

Product name	Mouse anti-Human blood group antigen B, clone HEB-29 (Monoclonal)
Host	Mouse
Applications	IHC-P, FUNC
Species reactivity	Human
Conjugate	-
Immunogen	Mixture of erythrocytes of group B and glycoprotein fraction isolated from saliva of secretors with blood group B.
Isotype	IgM
Clonality	Monoclonal
Clone number	HEB-29
Size	1.0 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with <0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human blood group antigen B, clone HEB-29 (Monoclonal)

Clone no. HEB-29

MONOSAN

Additional info

The antibody HEB-29 reacts with human blood group B. Agglutination test + B, A1B, A2B. The specificity of the antibody HEB-29 was confirmed by comparison of specificity and reactivity to standard reagent using >5 000 samples of blood.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES