

Mouse anti-Human blood group antigen A, clone HE-195 (Monoclonal)

Clone no. HE-195

MONOSAN

Product name	Mouse anti-Human blood group antigen A, clone HE-195 (Monoclonal)
Host	Mouse
Applications	FUNC
Species reactivity	Human
Conjugate	-
Immunogen	Mixture of erythrocytes of group A1 and glycoprotein fraction isolated from saliva of secretors with blood group A.
Isotype	IgM
Clonality	Monoclonal
Clone number	HE-195
Size	1.0 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with <0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human blood group antigen A, clone HE-195 (Monoclonal)

Clone no. HE-195

MONOSAN

Additional info

Human blood group A antigen belongs to a group of carbohydrate determinants carried on both glycolipids and glycoproteins; it is detected on erythrocytes and certain epithelial cells. The antibody HE-195 recognizes human blood group A including weak variants A3 , AX , A3B, AXB. The specificity of antibody HE-195 was confirmed by comparison of specificity and reactivity to standard reagent using >5.000 samples of blood. Agglutination test A1, A1B, A2B.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES